

To Whom It May Concern

Statement

Regarding the product Coronavirus (2019-nCoV) -Antigentest- produced by Beijing Hotgen Biotech Co., Ltd., the following statement are especially made:

The Coronavirus (2019-nCoV) -Antigentest- by our company detects N protein of SARS-CoV-2. While the mutation of the virus found in England (B.1.1.7 strain), South Africa (B.1.351 strain, **B.1.1.529 strain**), Brazil (P.1 strain) and India (B.1.617 strain, Delta variant B.1.617.2) are mainly on S protein. So our product is still effective to detect the antigen of novel coronavirus and the mutation will not affect the detection result of our product.

Beijing Hotgen Biotech Co., Ltd.

November 28, 2021

